



Fillmaster Digital

The Fillmaster 3600 is a **quick, error-proof way to dispense water**, dispensing exactly how much water is needed in a timely manner for reconstitutions. Pharmacists scan the NDC on the bottle for accurate and fast water dispensing.

Key Features

- Integrated reconstitution database
- NDC barcode reader (optional)
- Connects to FillPure Water Purification System
- Battery powered
- Push button operation
- Bright display

Specifications

- Dimensions: 18" x 4" x 3"
- Wall mountable
- Battery power: (4.5v)
- Dispensing quantity: 9mL to 999mL
- Accuracy: +/- 0.3mL at 9mL, +/- 3mL at 300mL