



FLAVORx Auto

Fully Automated Dispensing System

Every second counts in today's fast-paced, workflow-optimized pharmacy environment. The FLAVORx Auto has **completely eliminated the time it takes to custom-flavor liquid medications and other suspensions.**

By integrating completely into your work flow, the FLAVORx Auto allows your pharmacy to seamlessly reconstitute and add flavor into a medication in 15 seconds or less.

Why You'll Love it

1. Automated reconstitution
2. Automated flavoring
3. Barcode scanning
4. Wi-Fi connectivity
5. 5" touch screen display
6. Holds twelve 4oz flavor bottles
7. Connects to FillPure Water Purification System

Specifications

- Dimensions: 21.5" x 15.5" x 7"
- Weight: 22 lbs (fully loaded)

Inside the Box

- FLAVORx Auto dispenser
- ¼" white tubing: 50ft
- Power supply: 12v/6A power supply
- 12 FLAVORx bottles (120 mL) with custom caps
- SD card with NDC and barcode data
- Flavoring and reconstitution database
- FillPure Water Purification System
- FLAVORx marketing materials

